

PAK2 Rabbit mAb


ABZ2186

Product Overview


Name	PAK2 Rabbit mAb
Accession(Primary)	Q13177
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:20; IP-1:20-1:50

Target information (Q13177)

Gene Name	PAK2
About Gene Symbol	PAK2 (P21 (RAC1) Activated Kinase 2) is a serine/threonine protein kinase involved in cell cycle regulation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PAK2 in K562, rat Brain, C6, 3T3, Hela lysates using PAK2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using PAK2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.