

ERK5 Rabbit mAb


ABZ2187

Product Overview

Name	ERK5 Rabbit mAb
Accession(Primary)	Q13164
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 88 kDa; Observed MW: 115 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:20-1:50; IP-1:20-1:50


Target information (Q13164)

Gene Name	MAPK7
About Gene Symbol	



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of ERK5 in Hela lysates
using ERK5 antibody.



Western blot analysis of ERK5 in Hela lysates
using ERK5 antibody.