

FADD Rabbit mAb


ABZ2189

Product Overview

Name	FADD Rabbit mAb
Accession(Primary)	Q13158
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 23 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:20-1:100; ICC/IF-1:20-1:50; FC-1:20; IP-1:20-1:50

Target information (Q13158)

Gene Name	FADD
About Gene Symbol	



Western blot analysis of FADD in rat Brain, 3T3
lysates using FADD antibody.