


PAK1 Rabbit mAb

ABZ2190

Product Overview


Name	PAK1 Rabbit mAb
Accession(Primary)	Q13153
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20

Target information (Q13153)

Gene Name PAK1**About Gene Symbol**

Immunocytochemistry analysis of PAK1 (green) in Jurkat using PAK1 antibody, and DAPI(blue).





Western blot analysis of PAK1 in C6 lysates using PAK1 antibody.