

PAK1 Rabbit mAb


ABZ2191

Product Overview


Name	PAK1 Rabbit mAb
Accession(Primary)	Q13153
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 61 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:50-1:200; ICC/IF-1:50-1:100; FC-1:20; IP-1:20-1:50

Target information (Q13153)


Gene Name	PAK1
About Gene Symbol	




Immunocytochemistry analysis of PAK1 (green) in 3T3 using PAK1 antibody, and DAPI (blue).



Western blot analysis of PAK1 in 3T3 lysates using PAK1 antibody.



Western blot analysis of PAK1 in HeLa, CHO-K1, C6, Jurkat lysates using PAK1 antibody



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using PAK1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.