

TDP43 Rabbit mAb


ABZ2192

Product Overview


Name	TDP43 Rabbit mAb
Accession(Primary)	Q13148
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 45 kDa; Observed MW: 45 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:10

Target information (Q13148)


Gene Name	TARDBP
About Gene Symbol	TARDBP (TAR DNA Binding Protein) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of TDP 43 (green)
in 293t using TDP 43 antibody, and DAPI(blue).



Immunohistochemistry analysis of paraffin-
embedded Human brain using TDP 43
antibody. High-pressure and temperature Sodium
Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of TDP43 in K562, C6, 3T3, HeLa lysates using TDP43 antibody