

MTAP Rabbit mAb


ABZ2194

Product Overview


Name	MTAP Rabbit mAb
Accession(Primary)	Q13126
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 31 kDa; Observed MW: 31 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:50-1:100; FC-1:200-1:500

Target information (Q13126)

Gene Name	MTAP
About Gene Symbol	MTAP (Methylthioadenosine Phosphorylase) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of MTAP in 3T3 lysates using MTAP antibody.



Western blot analysis of MTAP in A549, MCF-7, U251, HL-60 lysates using MTAP antibody.