

IFT88 Rabbit mAb


ABZ2196

Product Overview


Name	IFT88 Rabbit mAb
Accession(Primary)	Q13099
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 93 kDa; Observed MW: 93 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:200

Target information (Q13099)

Gene Name	IFT88
About Gene Symbol	IFT88 (Intraflagellar Transport 88) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate



Western blot analysis of IFT88 in mouse testis lysates using IFT88 antibody.



Western blot analysis of IFT88 in Jurkat, C6, HeLa lysates using IFT88 antibody.