

MSP/MST1 Rabbit mAb


ABZ2198

Product Overview


Name	MSP/MST1 Rabbit mAb
Accession(Primary)	Q13043
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 56 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:200; ICC/IF-1:100-1:200; FC-1:20-1:50; IP-1:20-1:100

Target information (Q13043)

Gene Name	STK4
About Gene Symbol	STK4 (Serine/Threonine Kinase 4) is a serine/threonine protein kinase involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Immunocytochemistry analysis of MSP (green) in K562 using MSP antibody, and DAPI (blue).



Western blot analysis of MSP/MST1 in K562, rat Brain, C6, 3T3, HeLa lysates using MSP/MST1 antibody.