

BNIP3 Rabbit mAb


ABZ2199

Product Overview


Name	BNIP3 Rabbit mAb
Accession(Primary)	Q12983
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 22 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100

Target information (Q12983)


Gene Name	BNIP3
About Gene Symbol	BNIP3 (BCL2 Interacting Protein 3) is a structural cytoskeletal protein involved in apoptosis. It is an important research target in infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of BNIP3 in SH-SY5Y lysates using BNIP3 antibody.



Western blot analysis of BNIP3 in rat Brain lysates using BNIP3 antibody.



Immunocytochemistry analysis of BNIP3(green) in Hela using BNIP3 antibody, and DAPI(blue)