


PTP4A2 Rabbit mAb

ABZ2200

Product Overview

Name	PTP4A2 Rabbit mAb
Accession(Primary)	Q12974
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 19 kDa; Observed MW: 19 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q12974)

Gene Name PTP4A2**About Gene Symbol**

Western blot analysis of PTP4A2/PRL2 in HeLa, A549, HL-60 lysates using PTP4A2 antibody.

