


SAP97 Rabbit mAb

ABZ2201


Product Overview

Name	SAP97 Rabbit mAb
Accession(Primary)	Q12959
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 100 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q12959)

Gene Name DLG1**About Gene Symbol**

Western blot analysis of SAP97 in mouse brain, rat brain lysates using SAP97 antibody.



Western blot analysis of SAP97 in HeLa, A549, HL-60 lysates using SAP97 antibody.