

FOXC1 Rabbit mAb


ABZ2202

Product Overview

Name	FOXC1 Rabbit mAb
Accession(Primary)	Q12948
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 75 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:1000; FC-1:200-1:500; IP-1:20-1:50

Target information (Q12948)

Gene Name	FOXC1
About Gene Symbol	



Western blot analysis of FOXC1 in 3T3 lysates
using FOXC1 antibody.