

## ILF3 Rabbit mAb


ABZ2204

### Product Overview



|                             |                                                                           |
|-----------------------------|---------------------------------------------------------------------------|
| <b>Name</b>                 | ILF3 Rabbit mAb                                                           |
| <b>Accession(Primary)</b>   | Q12906                                                                    |
| <b>Host Species</b>         | Rabbit                                                                    |
| <b>Clonality</b>            | Monoclonal                                                                |
| <b>Isotype</b>              | IgG                                                                       |
| <b>Species Reactivity</b>   | Human                                                                     |
| <b>Calculated MW</b>        | Calculated MW: 95 kDa; Observed MW: 95 kDa                                |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:1000-1:5000;<br>ICC/IF-1:50-1:200; FC-1:20-1:50 |

### Target information (Q12906)


|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | ILF3 |
| <b>About Gene Symbol</b> |      |



Immunocytochemistry analysis of ILF3 (green) in Jurkat using ILF3 antibody, and DAPI (blue).



Western blot analysis of ILF3/NF90 in K562, 3T3, HeLa lysates using ILF3/NF90 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using ILF3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.