


SF3A3 Rabbit mAb

ABZ2206


Product Overview

Name	SF3A3 Rabbit mAb
Accession(Primary)	Q12874
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 59 kDa; Observed MW: 59 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:10-1:100

Target information (Q12874)

Gene Name SF3A3**About Gene Symbol**

Western blot analysis of SF3A3 in 3T3, Raw264.7, C6 lysates using SF3A3 antibody.



Western blot analysis of SF3A3 in HeLa, CHO-K1, C6, Jurkat lysates using SF3A3 antibody.