

Syntaxin 4 Rabbit mAb

ABZ2209


Product Overview

Name	Syntaxin 4 Rabbit mAb
Accession(Primary)	Q12846
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q12846)

Gene Name STX4

About Gene Symbol



Western blot analysis of Syntaxin 4 in C6, 3T3, HeLa lysates using Syntaxin 4 antibody.

