

SNF5/SMARCB1 Rabbit mAb


ABZ2210

Product Overview


Name	SNF5/SMARCB1 Rabbit mAb
Accession(Primary)	Q12824
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:50-1:100; FC-1:50-1:100; IP-1:20-1:50

Target information (Q12824)


Gene Name	SMARCB1
About Gene Symbol	



Western blot analysis of SNF5 in MCF-7, 3T3, HeLa lysates using SMARCB1 antibody.



Immunocytochemistry analysis of SNF5(green) in HeLa using SNF5 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using SNF5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.