

HspA14 Rabbit mAb


ABZ2212

Product Overview



Name	HspA14 Rabbit mAb
Accession(Primary)	Q0VDF9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200

Target information (Q0VDF9)


Gene Name	HSPA14
About Gene Symbol	



Immunocytochemistry analysis of HspA14 (green)
in Jurkat using HspA14 antibody, and DAPI(blue).



Western blot analysis of HSPA14 in 3T3 lysates
using HSPA14 antibody.



Western blot analysis of HSPA14 in HL-60, U251, U87-MG lysates using HSPA14 antibody