

LGALS3BP Rabbit mAb

ABZ2215


Product Overview

Name	LGALS3BP Rabbit mAb
Accession(Primary)	Q08380
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 65 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:500-1:1000; FC-1:100-1:200

Target information (Q08380)

Gene Name LGALS3BP

About Gene Symbol



Western blot analysis of LGALS3BP in HT-1080, HepG2, A549 lysates using LGALS3BP antibody.

