

Calcineurin A Rabbit mAb


ABZ2216

Product Overview


Name	Calcineurin A Rabbit mAb
Accession(Primary)	Q08209
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 59 kDa; Observed MW: 59 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:50-1:100

Target information (Q08209)

Gene Name	PPP3CA
About Gene Symbol	PPP3CA (Protein Phosphatase 3 Catalytic Subunit Alpha) is a protein phosphatase involved in cellular metabolism. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mec



Western blot analysis of Calcineurin A in rat brain, mouse lung lysates using Calcineurin A antibody.



Western blot analysis of Calcineurin A in HeLa, A549, U2OS lysates using Calcineurin A antibody.