

ACTN3 Rabbit mAb

ABZ2217


Product Overview

Name	ACTN3 Rabbit mAb
Accession(Primary)	Q08043
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 103 kDa; Observed MW: 103 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; FC-1:20-1:50

Target information (Q08043)


Gene Name	ACTN3
About Gene Symbol	

Rat skeletal muscle



Western blot analysis of ACTN3 in rat skeletal muscle lysates using ACTN3 antibody.

3T3



Western blot analysis of ACTN3 in 3T3 lysates using ACTN3 antibody.