

SF2 Rabbit mAb


ABZ2218

Product Overview


Name	SF2 Rabbit mAb
Accession(Primary)	Q07955
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 28 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:500; FC-1:50-1:100

Target information (Q07955)


Gene Name	SRSF1
About Gene Symbol	




Immunocytochemistry analysis of SF2 (green) in 293T using SF2 antibody, and DAPI (blue).



Western blot analysis of SF2 in 3T3 lysates using SF2 antibody.



Western blot analysis of SF2 in K562, C6 lysates using SF2 antibody



Immunohistochemistry analysis of paraffin-embedded Human tonsil using SF2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.