

MCL1 Rabbit mAb


ABZ2219

Product Overview


Name	MCL1 Rabbit mAb
Accession(Primary)	Q07820
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 37 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:500; FC-1:100-1:200

Target information (Q07820)

Gene Name	MCL1
About Gene Symbol	



Western blot analysis of MCL1 in C6 lysates using MCL1 antibody.



Immunocytochemistry analysis of MCL1 (green) in MCF-7 using MCL1 antibody, and DAPI (blue)