

## Bax Rabbit mAb


ABZ2221

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Bax Rabbit mAb                             |
| <b>Accession(Primary)</b>   | Q07812                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 21 kDa; Observed MW: 21 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IP-1:20-1:50              |

### Target information (Q07812)

|                          |                                                                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | BAX                                                                                                                                                                                                                                                                                                           |
| <b>About Gene Symbol</b> | <p>BAX (BCL2 Associated X, Apoptosis Regulator) is a protein kinase and signaling regulator involved in apoptosis. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.</p> |



Western blot analysis of Bax in C6, 3T3, HeLa  
lysates using Bax antibody.