

Bax Rabbit mAb


ABZ2222

Product Overview


Name	Bax Rabbit mAb
Accession(Primary)	Q07812
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:5000-1:30000; IHC-1:100-1:800; IP-1:20-1:100

Target information (Q07812)

Gene Name	BAX
About Gene Symbol	BAX (BCL2 Associated X, Apoptosis Regulator) is a protein kinase and signaling regulator involved in apoptosis. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Bax in C6, 3T3, HeLa lysates using Bax antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Bax antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.