

FMRP Rabbit mAb


ABZ2228

Product Overview


Name	FMRP Rabbit mAb
Accession(Primary)	Q06787
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 71 kDa; Observed MW: 70-80 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200

Target information (Q06787)


Gene Name	FMR1
About Gene Symbol	



Western blot analysis of FMRP in 293, U2OS, Jurkat, COS7 lysates using FMRP antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using FMRP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of FMRP (green) in
Hela using FMRP antibody, and DAPI(blue)