

RBPJK Rabbit mAb


ABZ2231

Product Overview

Name	RBPJK Rabbit mAb
Accession(Primary)	Q06330
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 56 kDa; Observed MW: 61 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100; FC-1:50-1:100

Target information (Q06330)

Gene Name	RBPJ
About Gene Symbol	



Western blot analysis of RBPJK in rat Brain, C6
lysates using RBPJK antibody.