


FPGS Rabbit mAb

ABZ2234


Product Overview

Name	FPGS Rabbit mAb
Accession(Primary)	Q05932
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 65 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:50-1:100

Target information (Q05932)

Gene Name FPGS**About Gene Symbol**

Western blot analysis of FPGS in rat lung lysates using FPGS antibody.



Western blot analysis of FPGS in K562, 3T3
lysates using FPGS antibody.