


NMDAR1 Rabbit mAb

ABZ2237

Product Overview

Name	NMDAR1 Rabbit mAb
Accession(Primary)	Q05586
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 105 kDa; Observed MW: 120 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50

Target information (Q05586)

Gene Name GRIN1**About Gene Symbol**

Western blot analysis of NMDAR1 in rat brain, Jurkat lysates using NMDAR1 antibody.

