



Ubiquitin Protein Ligase E3A Rabbit mAb

ABZ2240


Product Overview

Name	Ubiquitin Protein Ligase E3A Rabbit mAb
Accession(Primary)	Q05086
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 101 kDa; Observed MW: 101 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200; IP-1:20


Target information (Q05086)

Gene Name UBE3A**About Gene Symbol**

Immunocytochemistry analysis of Ubiquitin Protein Ligase E3A (green) in HeLa using Ubiquitin Protein Ligase E3A antibody, and DAPI(blue).



Western blot analysis of UBE3A in mouse brain, rat brain lysates using UBE3A antibody.



Western blot analysis of UBE3A in A549, HL-60, U2OS lysates using UBE3A antibody