


c-Rel Rabbit mAb

ABZ2241

Product Overview

Name	c-Rel Rabbit mAb
Accession(Primary)	Q04864
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 69 kDa; Observed MW: 78 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:10-1:100

Target information (Q04864)

Gene Name REL**About Gene Symbol**

Western blot analysis of c-Rel in HeLa, A549, HL-60 lysates using c-Rel antibody.

