


ACVR1 Rabbit mAb

ABZ2242


Product Overview

Name	ACVR1 Rabbit mAb
Accession(Primary)	Q04771
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 57 kDa; Observed MW: 57 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:10-1:100

Target information (Q04771)

Gene Name ACVR1**About Gene Symbol**

Western blot analysis of Activin Receptor Type IA in mouse heart, mouse brain lysates using Activin Receptor Type IA antibody.



Western blot analysis of Activin Receptor Type IA in HeLa lysates using Activin Receptor Type IA antibody.