


PKC theta Rabbit mAb

ABZ2243

Product Overview

Name	PKC theta Rabbit mAb
Accession(Primary)	Q04759
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 82 kDa; Observed MW: 79 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q04759)

Gene Name PRKCQ**About Gene Symbol**

Western blot analysis of PKC theta/PRKCQ in U2OS, C6 lysates using PKC theta antibody.

