

PSMA Rabbit mAb


ABZ2246

Product Overview

Name	PSMA Rabbit mAb
Accession(Primary)	Q04609
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 84 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; IP-1:10-1:100

Target information (Q04609)

Gene Name	FOLH1
About Gene Symbol	



Western blot analysis of PSMA in 3T3 lysates
using PSMA antibody.