

PSMA Rabbit mAb


ABZ2247

Product Overview


Name	PSMA Rabbit mAb
Accession(Primary)	Q04609
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 84 kDa; Observed MW: 110 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; FC-1:20; IP-1:50

Target information (Q04609)

Gene Name	FOLH1
About Gene Symbol	



Western blot analysis of PSMA in LnCaP lysates using PSMA antibody.



Immunohistochemistry analysis of paraffin-embedded Human prostate carcinoma using PSMA antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.