

PKR Rabbit mAb


ABZ2250

Product Overview


Name	PKR Rabbit mAb
Accession(Primary)	Q03963
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 74 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:100-1:200; FC-1:50-1:100; IP-1:20-1:50

Target information (Q03963)

Gene Name	Eif2Ak2
About Gene Symbol	



Immunocytochemistry analysis of PKR (green) in 3T3 using PKR antibody, and DAPI (blue).



Western blot analysis of PKR in C6, 3T3 lysates using PKR antibody.