

MVK Rabbit mAb


ABZ2251

Product Overview

Name	MVK Rabbit mAb
Accession(Primary)	Q03426
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:10-1:100

Target information (Q03426)

Gene Name	MVK
About Gene Symbol	MVK (Mevalonate Kinase) is a membrane receptor involved in translational regulation. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of MVK in HeLa, A549, HL-60, U251, U87-MG lysates using MVK antibody.