

CENPC Rabbit mAb


ABZ2253

Product Overview

Name	CENPC Rabbit mAb
Accession(Primary)	Q03188
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 107 kDa; Observed MW: 107 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:500; FC-1:200-1:500; IP-1:50-1:100


Target information (Q03188)

Gene Name	CENPC
About Gene Symbol	



Copyright© 2024 Zen-bioscience Limited.

Western blot analysis of CENPC in Jurkat and K562 lysates using CENPC antibody.



Western blot analysis of CENPC in K562 lysates using CENPC antibody.