

Caveolin-1 Rabbit mAb


ABZ2254

Product Overview


Name	Caveolin-1 Rabbit mAb
Accession(Primary)	Q03135
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 20 kDa; Observed MW: 20 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:500-1:2000; ICC/IF-1:100-1:200; FC-1:50-1:100; IP-1:20-1:50

Target information (Q03135)


Gene Name	CAV1
About Gene Symbol	CAV1 (Caveolin 1) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Caveolin1 in HeLa lysates using Caveolin 1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Caveolin1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of Caveolin1
(green) in Hela using Caveolin1 antibody, and
DAPI(blue)