

FKBP52 Rabbit mAb


ABZ2255

Product Overview


Name	FKBP52 Rabbit mAb
Accession(Primary)	Q02790
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100

Target information (Q02790)


Gene Name	FKBP4
About Gene Symbol	



Western blot analysis of FKBP52 in Raw264.7 lysates using FKBP52 antibody.



Western blot analysis of FKBP52 in Jurkat, C6, CHO-K1, HeLa lysates using FKBP52 antibody.



Immunocytochemistry analysis of FKBP52 (green)
in hela using FKBP52 antibody, and DAPI(blue)