


IRF8 Rabbit mAb

ABZ2258

Product Overview

Name	IRF8 Rabbit mAb
Accession(Primary)	Q02556
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 48 kDa; Observed MW: 48 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IP-1:20

Target information (Q02556)

Gene Name IRF8**About Gene Symbol**

Western blot analysis of IRF8 in THP-1 lysates using IRF8 antibody.

