

## MKLP1 Rabbit mAb


ABZ2259

### Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | MKLP1 Rabbit mAb                             |
| <b>Accession(Primary)</b>   | Q02241                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                            |
| <b>Calculated MW</b>        | Calculated MW: 110 kDa; Observed MW: 110 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; ICC/IF-1:200-1:500          |

### Target information (Q02241)

|                          |                                                                                                                                                                                                                                                                                              |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | KIF23                                                                                                                                                                                                                                                                                        |
| <b>About Gene Symbol</b> | KIF23 (Kinesin Family Member 23) is a kinesin motor protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of MKLP1 in HeLa lysates  
using MKLP1 antibody.