


CENPE Rabbit mAb

ABZ2260

Product Overview

Name	CENPE Rabbit mAb
Accession(Primary)	Q02224
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 316 kDa; Observed MW: 316 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q02224)

Gene Name CENPE**About Gene Symbol**

Western blot analysis of CENPE in Jurkat lysates using CENPE antibody.

