

MEF2A Rabbit mAb

ABZ2261


Product Overview

Name	MEF2A Rabbit mAb
Accession(Primary)	Q02078
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:2000-1:5000; ICC/IF-1:100-1:200


Target information (Q02078)

Gene Name MEF2A

About Gene Symbol



Western blot analysis of MEF2A in K562, rat Brain, C6, 3T3, HeLa lysates using MEF2A antibody.



Immunocytochemistry analysis of MEF2A (green)
in MCF-7 using MEF2A antibody, and DAPI (blue)