

EWSR1 Rabbit mAb


ABZ2263

Product Overview


Name	EWSR1 Rabbit mAb
Accession(Primary)	Q01844
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 68 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:200-1:1000; FC-1:50-1:100

Target information (Q01844)


Gene Name	EWSR1
About Gene Symbol	EWSR1 (EWS RNA Binding Protein 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of EWSR1 in rat brain lysates using EWSR1 antibody.



Western blot analysis of EWSR1 in MCF-7, 3T3, HeLa lysates using EWSR1 antibody.



Immunocytochemistry analysis of EWSR1(green)
in Hela using EWSR1 antibody, and DAPI(blue)