

PFKP Rabbit mAb


ABZ2265

Product Overview


Name	PFKP Rabbit mAb
Accession(Primary)	Q01813
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 86 kDa; Observed MW: 86 kDa
Recommended Dilution	WB-1:1000-1:10000; ICC/IF-1:100-1:200; FC-1:100-1:200; IP-1:50-1:100

Target information (Q01813)

Gene Name	PFKP
About Gene Symbol	PFKP (Phosphofructokinase, Platelet) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PFKP in C6, 3T3, HeLa lysates using PFKP antibody.



Immunocytochemistry analysis of PFKP(green) in HeLa using PFKP antibody, and DAPI(blue)