

Fragilis Rabbit mAb


ABZ2266

Product Overview


Name	Fragilis Rabbit mAb
Accession(Primary)	Q01628
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; ICC/IF-1:50-1:100

Target information (Q01628)


Gene Name	IFITM3
About Gene Symbol	



Western blot analysis of Fragilis in K562, HeLa lysates using Fragilis antibody.



Immunocytochemistry analysis of Fragilis(green) in HeLa using Fragilis antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using Fragilis antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.