

FLI1 Rabbit mAb


ABZ2267

Product Overview


Name	FLI1 Rabbit mAb
Accession(Primary)	Q01543
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 51 kDa; Observed MW: 51 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:200-1:500

Target information (Q01543)

Gene Name	FLI1
About Gene Symbol	FLI1 (Fli-1 Proto-Oncogene, ETS Transcription Factor) is a proto-oncogene involved in gene transcription and regulation. It is an important research target in cancer research and infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of FLI1 in Jurkat, C6 lysates using FLI1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using FLI1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.