

SET Rabbit mAb


ABZ2270

Product Overview


Name	SET Rabbit mAb
Accession(Primary)	Q01105
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:500-1:1000; FC-1:10-1:100; IP-1:20-1:50

Target information (Q01105)


Gene Name	SET
About Gene Symbol	



Western blot analysis of SET in K562, HeLa lysates using SET antibody.



Immunocytochemistry analysis of SET(green) in HeLa using SET antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using SET antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.