

E2F1 Rabbit mAb


ABZ2271

Product Overview


Name	E2F1 Rabbit mAb
Accession(Primary)	Q01094
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:10000; IHC-1:50-1:500; IHC-1:50-1:100; ICC/IF-1:50-1:100

Target information (Q01094)


Gene Name	E2F1
About Gene Symbol	E2F1 (E2F Transcription Factor 1) is a transcription factor involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.




Western blot analysis of E2F1 in MOLT4 lysates using E2F1 antibody.



Western blot analysis of E2F1 in C6, 3T3 lysates using E2F1 antibody.



Immunocytochemistry analysis of E2F1(green) in Hela using E2F1 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human tonsil using E2F1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.